


# PAS3000 Porosity Analysis System by Porosity Solutions

- ★ Control the hydrogen porosity in your melt
- ★ Scientific method of evaluating RPT Samples (gas slug)
- ★ No need to saw and polish RPT sample - Safer, cheaper, more accurate! -
- ★ Enclosed for foundry floor operation
- ★ Automatic Calculation of Specific Gravity
- ★ Wide - easy access tray
- ★ Results in just a few seconds!



Precision Calibration weights included!

# PAS3000 Porosity Analysis System by Porosity Solutions

*Designed and Built for the Shop Floor! – Rugged, Easy to Use, Accurate*

**LED illuminated interior  
for easy viewing**

**Fast 3 Step Controls**

1. Tare
2. Weigh in Air
3. Weigh in Water

**The Result is  
automatically displayed**

**Rugged wide load tray  
for easy sample load  
and unload**

**Ethernet port for  
computer or network  
communication**

**Industrial Ohaus scale  
that communicates with  
the control panel**






# Easy as 1 – 2 – 3



Weigh in Air



Weigh in Water



Result

Porosity Solutions - Contact Brad Hohenstein  
blh@porositysolutions.com